

JRKL Antibody

Catalog No: #47677

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	JRKL Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antigen affinity purification
Applications	WB
Species Reactivity	Hu, Ms
Specificity	The antibody detects endogenous levels of total JRKL protein.
Immunogen Type	Peptide
Immunogen Description	Synthetic peptide of human JRKL
Target Name	JRKL
Other Names	HHMJG
Accession No.	Swiss-Prot#:Q9Y4A0NCBI Gene ID:8690Gene Accssion:NP_003763
Calculated MW	60 kDa
Concentration	1
Formulation	Rabbit IgG in pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol.
Storage	Store at -20°C

Application Details

WB dilution:1:1000-5000

Background

The function of this gene has not yet been defined, however, the encoded protein shares similarity with the human (41% identical) and mouse (34% identical) jerky gene products. This protein may act as a nuclear regulatory protein. Alternatively spliced transcript variants have been found for this gene.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.