

DYKDDDDK Tag(FLAG) HRP Antibody

Catalog No: #48045



Package Size: #48045-1 50ul #48045-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	DYKDDDDK Tag(FLAG) HRP Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Immunogen affinity purified
Applications	WB
Immunogen Description	peptide
Conjugates	HRP-conjugated
Other Names	DDDDK epitope tag antibody DDDK antibody ddk antibody DYKDDDDK antibody DYKDDDDK epitope tag antibody DYKDDDDK tag antibody ECS epitope tag antibody ECS tag antibody Enterokinase Cleavage Site epitope tag antibody Enterokinase Cleavage Site tag antibody FLAG antibody FLAG tag antibody
Accession No.	Swiss-Prot#:
Formulation	1*TBS (pH7.4), 1%BSA, 40%Glycerol. Preservative: 0.05% Sodium Azide.
Storage	Store at -20°C


Application Details

WB: 1:5,000

Images



Western blot analysis on DYKDDDDK Tag recombinant protein (80kDa).



Background

The DYKDDDDK peptide has been used extensively as a general epitope tag in expression vectors. Due to its small size, it is unlikely to affect the tagged protein's biochemical properties. This antibody recognizes the DYKDDDDK peptide fused to either the amino- or carboxy-terminus of targeted proteins.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.