

V5-Tag HRP Antibody

Catalog No: #48053

Package Size: #48053-1 50ul #48053-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	V5-Tag HRP Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Immunogen affinity purified
Applications	WB
Immunogen Description	peptide
Conjugates	HRP-conjugated
Other Names	GKPIPPLLGLDST epitope tag antibody GKPIPPLLGLDST tag antibody Protein Rev antibody Regulator of expression of viral proteins antibody rev antibody V5 epitope tag antibody
Accession No.	Swiss-Prot#:
Calculated MW	50 kDa
Formulation	1*TBS (pH7.4), 1%BSA, 40%Glycerol. Preservative: 0.05% Sodium Azide.
Storage	Store at -20°C

Application Details

WB: 1:10,000o½ 1:100,000

Images



Western blot analysis on V5-tagged recombinant protein (50kDa) using anti-V5 polyclonal antibody.

Background

The V5 epitope tag is derived from a small epitope (Pk) present on the P and V proteins of the paramyxovirus of simian virus 5 (SV5). The anti-V5 antibody detects recombinant proteins containing the 14 amino acid GKPIPPLLGLDST sequence.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.