

# Transmembrane protein 200A Antibody

Catalog No: #48081

Package Size: #48081-1 50ul #48081-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Transmembrane protein 200A Antibody
Host Species	Mouse
Clone No.	4-B3
Purification	Protein affinity purified
Applications	WB, ICC, IHC
Species Reactivity	Hu, Ms
Immunogen Description	Recombinant protein within human Transmembrane protein 200A 50-400 aa.
Other Names	KIAA1913 antibody T200A_HUMAN antibody Tmem200a antibody Transmembrane protein 200A antibody TTMC antibody Two transmembrane C antibody
Accession No.	Swiss-Prot#:Q86VY9
Calculated MW	47 kDa
Formulation	1*TBS (pH7.4), 1%BSA, 40%Glycerol. Preservative: 0.05% Sodium Azide.
Storage	Store at -20°C

## Application Details

WB: 1:500

IHC: 1:50

ICC: 1:100

## Background

Transmembrane protein 200A(TM200A) is multi-pass membrane protein, which belongs to the TM200 family. It is expressed in cerebellum.

## References

- Hui A B Y et al. Identification of a novel homozygous deletion region at 6q23.1 in medulloblastomas using high-resolution array comparative genomic hybridization analysis. Clin Cancer Res 11:4707-4716 (2005).
- Carninci P et al. The transcriptional landscape of the mammalian genome. Science 309:1559-1563 (2005).

Note: This product is for in vitro research use only and is not intended for use in humans or animals.