

SMCC Antibody

Catalog No: #48122



Package Size: #48122-1 50ul #48122-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	SMCC Antibody
Host Species	Mouse
Clone No.	6-F6
Purification	ProA affinity purified
Applications	ELISA
Species Reactivity	Species independent
Immunogen Description	Recombinant protein fused with SMCC.
Other Names	22 kDa actin binding protein antibody 22 kDa actin-binding protein antibody Human 22kDa smooth muscle protein antibody Protein WS3-10 antibody SM22 alpha antibody SM22 antibody SM22-alpha antibody SMCC antibody Smooth muscle protein 22-alpha antibody TAGL_HUMAN antibody TAGLN antibody TAGLN1 antibody Transgelin antibody Transgelin variant 2 antibody WS3 10 antibody
Accession No.	Swiss-Prot#:
Formulation	1*PBS (pH7.4), 0.2% BSA, 40% Glycerol. Preservative: 0.05% Sodium Azide.
Storage	Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Application Details

ELISAoO 1:50,000-1:100,000

References

1. Mattson, G., et al. (1993). A practical approach to cross-linking. *Molecular Biology Reports* 17:167-83.
2. Brinkley, M.A. (1992). A survey of methods for preparing protein conjugates with dyes, haptens and cross-linking reagents. *Bioconjugate Chem* 3:2-13.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.