

VSV-g Tag Antibody

Catalog No: #48191



Package Size: #48191-1 50ul #48191-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	VSV-g Tag Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Immunogen affinity purified
Applications	WB
Immunogen Description	peptide
Other Names	Vesicular stomatitis virus glycoprotein tag antibody VSV epitope tag antibody VSV tag antibody VSV-G epitope tag antibody VSVG antibody vsvg tag antibody
Accession No.	Swiss-Prot#:
Formulation	1*TBS (pH7.4), 1%BSA, 40%Glycerol. Preservative: 0.05% Sodium Azide.
Storage	Store at -20°C

Application Details

WB: 1:5,000-1:10,000

Images



Western blot analysis on VSV-g tagged fusion protein.

Background

VSV-G, a vesicular stomatitis virus G (VSV-G) protein fragment. VSV-G protein is commonly used in biomedical research to pseudotype retroviral and lentiviral vectors, conveying the ability to transduce a broad range of mammalian cell types with genes of interest. The VSV-G protein has also been used in cytological studies of trafficking in the endomembrane system. This antibody recognizes the VSV-g-tag fused to either the amino- or carboxy-terminus of targeted proteins.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.