

## SREBP-1 Polyclonal Antibody

Catalog No: #41458

Package Size: #41458-1 50ul #41458-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	SREBP-1 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC ELISA IF
Species Reactivity	Hu, Ms, Rt
Specificity	SREBP-1 Polyclonal Antibody detects endogenous levels of SREBP-1 protein.
Immunogen Description	Synthesized peptide derived from human SREBP-1 around the non-phosphorylation site of S439.
Target Name	SREBP-1
Other Names	SREBF1; BHLHD1; SREBP1; Sterol regulatory element-binding protein 1; SREBP-1; Class D basic helix-loop-helix protein 1; bHLHD1; Sterol regulatory element-binding transcription factor 1
Accession No.	Swiss-Prot: P36956NCBI Gene ID: 6720
SDS-PAGE MW	122kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

## Application Details

Western Blot: 1/500 - 1/2000.


Immunohistochemistry: 1/100 - 1/300.

ELISA: 1/20000.

IF 1:50-200

Not yet tested in other applications.

## Images



Western Blot analysis of hela cells using SREBP-1 Polyclonal Antibody

---

Note: This product is for in vitro research use only and is not intended for use in humans or animals.