

VASP Polyclonal Antibody

Catalog No: #41529



Package Size: #41529-1 50ul #41529-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	VASP Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC IF ELISA
Species Reactivity	Hu Ms Rt Mk
Specificity	VASP Polyclonal Antibody detects endogenous levels of VASP protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from human VASP around the non-phosphorylation site of S157.
Target Name	VASP
Other Names	VASP; Vasodilator-stimulated phosphoprotein; VASP
Accession No.	Swiss-Prot: P50552NCBI Gene ID: 7408
SDS-PAGE MW	46+50kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

Western Blot: 1/500 - 1/2000.


Immunohistochemistry: 1/100 - 1/300.

Immunofluorescence: 1/200 - 1/1000.

ELISA: 1/20000.

Not yet tested in other applications.

Images



Western Blot analysis of COS7 cells using VASP Polyclonal Antibody

Published Papers

el at., Identification of Novel Biomarkers With Diagnostic Value and Immune Infiltration in Burn Injury In Front Genet on 2022 Mar 22 by

Sitong Zhou, Kangchun Wang,et al..PMID:35391790, , (2022)

[PMID:35391790](#)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.