

## VCP Polyclonal Antibody

Catalog No: #41530



Package Size: #41530-1 50ul #41530-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	VCP Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC ELISA
Species Reactivity	Hu Ms Rt
Specificity	VCP Polyclonal Antibody detects endogenous levels of VCP protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from human VCP around the non-phosphorylation site of S352.
Target Name	VCP
Other Names	VCP; Transitional endoplasmic reticulum ATPase; TER ATPase; 15S Mg(2+)-ATPase p97 subunit; Valosin-containing protein; VCP
Accession No.	Swiss-Prot: P55072NCBI Gene ID: 7415
SDS-PAGE MW	85kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

## Application Details


Western Blot: 1/500 - 1/2000.

Immunohistochemistry: 1/100 - 1/300.

ELISA: 1/40000.

Not yet tested in other applications.

## Images



Western Blot analysis of K562 cells using VCP Polyclonal Antibody

---

Note: This product is for in vitro research use only and is not intended for use in humans or animals.