

XRCC1 Polyclonal Antibody

Catalog No: #41538

Package Size: #41538-1 50ul #41538-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	XRCC1 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC ELISA
Species Reactivity	Hu Ms Rt
Specificity	XRCC1 Polyclonal Antibody detects endogenous levels of XRCC1 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the C-terminal region of human XRCC1.
Target Name	XRCC1
Other Names	XRCC1; DNA repair protein XRCC1; X-ray repair cross-complementing protein 1
Accession No.	Swiss-Prot: P18887NCBI Gene ID: 7515
SDS-PAGE MW	70kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details


Western Blot: 1/500 - 1/2000.

Immunohistochemistry: 1/100 - 1/300.

ELISA: 1/10000.

Not yet tested in other applications.

Images



Western Blot analysis of NIH-3T3 Jurkat cells using XRCC1 Polyclonal Antibody

Note: This product is for in vitro research use only and is not intended for use in humans or animals.