

LXR α Polyclonal Antibody

Catalog No: #41558



Package Size: #41558-1 50ul #41558-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	LXR α Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB ELISA IHC IF
Species Reactivity	Hu Ms Rt
Specificity	LXR α Polyclonal Antibody detects endogenous levels of LXR α protein.
Immunogen Description	Synthesized peptide derived from the Internal region of human LXR α .
Target Name	LXR α
Other Names	NR1H3; LXRA; Oxysterols receptor LXR-alpha; Liver X receptor alpha; Nuclear receptor subfamily 1 group H member 3
Accession No.	Swiss-Prot: Q13133NCBI Gene ID: 10062
SDS-PAGE MW	50kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

Western Blot: 1/500 - 1/2000.


ELISA: 1/20000.


Not yet tested in other applications.

IHC-p: 1:100-300 .

IF 1:50-200

Images

Western Blot analysis of HeLa cells using LXR α Polyclonal Antibody



Immunohistochemical analysis of paraffin-embedded mouse-brain, antibody was diluted at 1:100

Note: This product is for in vitro research use only and is not intended for use in humans or animals.