

IL-5R α Polyclonal Antibody

Catalog No: #41622

Package Size: #41622-1 50ul #41622-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	IL-5R α Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB;IHC;IF;ELISA
Species Reactivity	Human;Rat;Mouse;
Specificity	IL-5R α Polyclonal Antibody detects endogenous levels of IL-5R α protein.
Immunogen Description	Synthesized peptide derived from the N-terminal region of human IL-5R α .
Target Name	IL-5R α
Other Names	IL5RA; IL5R; Interleukin-5 receptor subunit alpha; IL-5 receptor subunit alpha; IL-5R subunit alpha; IL-5R-alpha; IL-5RA; CDw125; CD125
Accession No.	Swiss-Prot: Q01344NCBI Gene ID: 3568
SDS-PAGE MW	47kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details


WB 1:500 - 1:2000.

IHC: 1:100-300

ELISA: 1:20000..

IF 1:50-200

Images

Western Blot analysis of K562 cells using IL-5R α Polyclonal Antibody

Note: This product is for in vitro research use only and is not intended for use in humans or animals.