

CD73 Polyclonal Antibody

Catalog No: #41661



Package Size: #41661-1 50ul #41661-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	CD73 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB ELISA IHC IF
Species Reactivity	Human, Mouse, Rat
Specificity	CD73 Polyclonal Antibody detects endogenous levels of CD73 protein.
Immunogen Description	Synthesized peptide derived from the Internal region of human CD73.
Target Name	CD73
Other Names	NT5E; NT5; NTE; 5'-nucleotidase; 5'-NT; Ecto-5'-nucleotidase; CD73
Accession No.	Swiss-Prot: P21589NCBI Gene ID: 4907
SDS-PAGE MW	68kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

Western Blot: 1/500 - 1/2000.


ELISA: 1/20000.

Not yet tested in other applications.


IHC: 1:100-300

IF 1:50-200

Images



Western Blot analysis of MCF7 L929 K562 cells using CD73 Polyclonal Antibody



Immunohistochemical analysis of paraffin-embedded human-skin, antibody was diluted at 1:100

Note: This product is for in vitro research use only and is not intended for use in humans or animals.