

## Nectin 3 Polyclonal Antibody

Catalog No: #41668



Package Size: #41668-1 50ul #41668-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Nectin 3 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB ELISA
Species Reactivity	Hu Ms Rt
Specificity	Nectin 3 Polyclonal Antibody detects endogenous levels of Nectin 3 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the Internal region of human Nectin 3.
Target Name	Nectin 3
Other Names	PVRL3; PRR3; Poliovirus receptor-related protein 3; CDw113; Nectin-3; CD113
Accession No.	Swiss-Prot: Q9NQS3NCBI Gene ID: 25945
SDS-PAGE MW	60kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


## Application Details

Western Blot: 1/500 - 1/2000.

ELISA: 1/20000.

Not yet tested in other applications.

## Images



Western Blot analysis of HeLa cells using Nectin 3 Polyclonal Antibody

## Published Papers

---

el et al., CD155 Cooperates with PD-1/PD-L1 to Promote Proliferation of Esophageal Squamous Cancer Cells via PI3K/Akt and MAPK Signaling Pathways. In *Cancers (Basel)* on 2022 Nov 15 by Xiyang Tan, Jie Yang, et al.. PMID:36428703, , (2022)

[PMID:36428703](#)

---

---

Note: This product is for in vitro research use only and is not intended for use in humans or animals.