

## IGFBP6 Polyclonal Antibody

Catalog No: #41750



Package Size: #41750-1 50ul #41750-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	IGFBP6 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB ELISA
Species Reactivity	Human;Rat;Mouse
Specificity	IGFBP6 Polyclonal Antibody detects endogenous levels of IGFBP6 protein.
Immunogen Description	Synthesized peptide derived from the Internal region of human IGFBP6.
Target Name	IGFBP6
Other Names	IGFBP6; IBP6; Insulin-like growth factor-binding protein 6; IBP-6; IGF-binding protein 6; IGFBP-6
Accession No.	Swiss-Prot: P24592NCBI Gene ID: 3489
SDS-PAGE MW	25kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


## Application Details

Western Blot: 1/500 - 1/2000.

ELISA: 1/20000.

Not yet tested in other applications.

## Images



Western Blot analysis of K562 cells using IGFBP6 Polyclonal Antibody.

---

Note: This product is for in vitro research use only and is not intended for use in humans or animals.