

LIR-7 Polyclonal Antibody

Catalog No: #41838



Package Size: #41838-1 50ul #41838-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	LIR-7 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB ELISA
Species Reactivity	Hu
Specificity	LIR-7 Polyclonal Antibody detects endogenous levels of LIR-7 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the Internal region of human LIR-7.
Target Name	LIR-7
Other Names	LILRA2; ILT1; LIR7; Leukocyte immunoglobulin-like receptor subfamily A member 2; CD85 antigen-like family member H; Immunoglobulin-like transcript 1; ILT-1; Leukocyte immunoglobulin-like receptor 7; LIR-7; CD85h
Accession No.	Swiss-Prot: Q8N149NCBI Gene ID: 11027
SDS-PAGE MW	53kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


Application Details

Western Blot: 1/500 - 1/2000.

ELISA: 1/20000.

Not yet tested in other applications.

Images



Western Blot analysis of HepG2 cells using LIR-7 Polyclonal Antibody.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.