

AChR α 1 Polyclonal Antibody

Catalog No: #41921

Package Size: #41921-1 50ul #41921-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	AChR α 1 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB ELISA IF
Species Reactivity	Hu Ms Rt
Specificity	AChR α 1 Polyclonal Antibody detects endogenous levels of AChR α 1 protein.
Immunogen Description	Synthesized peptide derived from the Internal region of human AChR α 1.
Target Name	AChR α 1
Other Names	CHRNA1; ACHRA; CHNRA; Acetylcholine receptor subunit alpha
Accession No.	Swiss-Prot: P02708NCBI Gene ID: 1134
SDS-PAGE MW	55kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details


Western Blot: 1/500 - 1/2000.

IF 1:100-300

ELISA: 1/10000.


Not yet tested in other applications.

Images



Western Blot analysis of NIH-3T3 cells using AChR α 1 Polyclonal Antibody

Immunofluorescence analysis of A549.



Note: This product is for in vitro research use only and is not intended for use in humans or animals.