

## AHSG Monoclonal Antibody

Catalog No: #42014

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

Product Name	AHSG Monoclonal Antibody
Host Species	Mouse
Clonality	Monoclonal
Purification	protein G purified
Applications	IHC
Species Reactivity	Hu
Specificity	specific for Human AHSG denatured and native forms
Immunogen Description	Recombinant Human AHSG
Target Name	AHSG
Other Names	Alpha-2-Z-globulin, Ba-alpha-2-glycoprotein, Fetuin-A, AHSG
Accession No.	Swiss-Prot#: P02765
Concentration	1.0mg/mL
Formulation	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Storage	Store at -20°C


## Application Details

Immunohistochemistry: 1:20 - 1:200

## Images



Immunohistochemical analysis of paraffin-embedded human liver cancer using #42014 at dilution of 1:200.



Immunohistochemical analysis of paraffin-embedded human breast cancer using #42014 at dilution of 1:200.

## Background

---

Promotes endocytosis, possesses opsonic properties and influences the mineral phase of bone. Shows affinity for calcium and barium ions.

---

Note: This product is for in vitro research use only and is not intended for use in humans or animals.