

## PGC Monoclonal Antibody

Catalog No: #42044



Package Size: #42044-1 50ul #42044-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

Product Name	PGC Monoclonal Antibody
Host Species	Mouse
Clonality	Monoclonal
Purification	protein G purified
Applications	IHC
Species Reactivity	Hu
Specificity	specific for PGC native and denatured forms
Immunogen Type	protein
Immunogen Description	Recombinant Human PGC
Target Name	PGC
Other Names	Pepsinogen C , PGC
Accession No.	Swiss-Prot#: P20142
Concentration	1.0mg/mL
Formulation	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Storage	Store at -20°C

## Application Details

Immunohistochemistry: 1:20 - 1:200

## Images



Immunohistochemical analysis of paraffin-embedded human gastric cancer tissue using #42044 at dilution of 1:200.

## Background

Hydrolyzes a variety of proteins.

---

Note: This product is for in vitro research use only and is not intended for use in humans or animals.