

MALT1 antibody

Catalog No: #22761

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description


Product Name	MALT1 antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Purified by antigen-affinity chromatography.
Applications	WB
Species Reactivity	Hu Ms Rt
Immunogen Description	Recombinant protein of human MALT1
Target Name	MALT1
Accession No.	NCBI Gene ID: 10892NCBI mRNA#: NM_173844NCBI Protein#: NP_776216
Concentration	1mg/ml
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.

Application Details

Predicted MW: 91kd

Western blotting: 1:500-1:2000

Images



Western blot analysis of extracts of hela cell lines, using MALT1 antibody. The lane on the left is treated with the antigen-specific peptide.

Background

This gene has been found to be recurrently rearranged in chromosomal translocation with two other genes - baculoviral IAP repeat-containing protein 3 (also known as apoptosis inhibitor 2) and immunoglobulin heavy chain locus - in mucosa-associated lymphoid tissue lymphomas. The protein encoded by this gene may play a role in NF-kappaB activation. Two alternatively spliced transcript variants encoding different isoforms have been described for this gene. [provided by RefSeq]

Note: This product is for in vitro research use only and is not intended for use in humans or animals.