

Vesicular, overexpressed in cancer, prosurvival protein 1 Polyclonal Antibody

Catalog No: #42397

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	Vesicular, overexpressed in cancer, prosurvival protein 1 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Caprylic Acid Ammonium Sulfate Precipitation purified
Applications	WB
Species Reactivity	Hu
Specificity	The antibody detects endogenous level of total Vesicular, overexpressed in cancer, prosurvival protein 1 polyclonal antibody.
Immunogen Type	protein
Immunogen Description	Recombinant human Vesicular, overexpressed in cancer, prosurvival protein 1 protein
Target Name	Vesicular, overexpressed in cancer, prosurvival protein 1
Other Names	VOPP1
Accession No.	Swiss-Prot#: Q96AW1
Calculated MW	19kd
Formulation	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Storage	Store at -20°C

Application Details

Western blotting: □1:500 - 1:1000

Images



All lanes : Vesicular, overexpressed in cancer, prosurvival protein 1 antibody at at 2ug/ml
 Lane 1 : EC109 whole cell lysate
 Lane 2 : 293T whole cell lysate
 Secondary Goat polyclonal to Rabbit IgG at 1/15000 dilution
 Predicted band size : 19 kDa
 Observed band size : 50 kDa

Background

Increases the transcriptional activity of NFKB1 by facilitating its nuclear translocation, DNA-binding and associated apoptotic response, when overexpressed.

References

[1] "Towards a catalog of human genes and proteins: sequencing and analysis of 500 novel complete protein coding human cDNAs." Wiemann S., Weil B., Wellenreuther R., Gassenhuber J., Glassl S., Ansorge W., Boecher M., Bloecker H., Bauersachs S., Blum H., La

Note: This product is for in vitro research use only and is not intended for use in humans or animals.