

Nod1 antibody

Catalog No: #22788



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Nod1 antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Purified by antigen-affinity chromatography.
Applications	WB IF
Species Reactivity	Hu
Immunogen Type	Recombinant protein
Immunogen Description	Recombinant protein fragment contain a sequence corresponding to a region within amino acids 698 and 928 of Human NOD1
Target Name	Nod1
Accession No.	NCBI Gene ID: 10392NCBI mRNA#: NM_006092NCBI Protein#: NP_006083
Concentration	1mg/ml
Formulation	Supplied in 0.1M Tris-buffered saline with 10% Glycerol (pH7.0). 0.01% Thimerosal was added as a preservative.
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.


Application Details

Predicted MW: 108kd

Western blotting: 1:500-1:3000

Immunofluorescence: 1:100-1:200

Images




Sample (30 ug of whole cell lysate)

A: H1299


7.5% SDS PAGE

Primary antibody diluted at 1: 1000

Immunofluorescence analysis of paraformaldehyde-fixed HeLa, using NOD1 antibody at 1: 200 dilution.



Merged with DNA probe



Note: This product is for in vitro research use only and is not intended for use in humans or animals.