

Mouse IgG

Catalog No: #IG101



Package Size: #IG101-1 100mg #IG101-2 200mg #IG101-3 1g

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Mouse IgG
Purification	Affinity chromatography or caprylic acid-ammonium sulfate,β 95%
Conjugates	None
Concentration	mg/ml
Storage	-20°C

Product Description

Source: Normal, healthy, unimmunized mouse serum

Appearance: Liquid

Filtration: 0.2um

Buffer :10mM Phosphate buffered salinePBSpH7.40.

Preservative: None

Stability: Stable for 24 months

Precautions: This product can be kept at 2-8Λ C for maximum 30 days. Avoid repeated freezing and thawing. If multiple usages are needed, it is recommended to be stored separately in smaller volumes.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.