

# Goat IgG

Catalog No: #IG501



Package Size: #IG501-1 100mg #IG501-2 200mg #IG501-3 1g

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Goat IgG
Purification	Affinity chromatography or caprylic acid-ammonium sulfate, $\beta$ 95%
Conjugates	None
Concentration	$\beta$ 20 mg/ml
Storage	-20°C

## Product Description

Source: Normal, healthy, unimmunized goat serum

Appearance: Liquid

Filtration: 0.2um

Buffer :10mM sodium Phosphate150mM sodium chloridePBSpH7.4

Preservative: None

Stability: Stable for 24 months

Precautions: This product can be kept at 2-8 $\Lambda$  C for maximum 30 days. Avoid repeated freezing and thawing. If multiple usages are needed, it is recommended to be stored separately in smaller volumes.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.