

CD137(human)(rec.)(His)

Catalog No: #80305



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	CD137(human)(rec.)(His)
Host Species	E. coli
Purification	>95%
Species Reactivity	Hu
Other Names	T-cell antigen ILA; CD137 antigen; CDw137; ILA; 4-1BB; MGC2172; 4-1BBR; TNFRSF9.
Calculated MW	19 KDa
Target Sequence	Human CD137[TNFRSF9/4-1BB] (AAH06196.1)(Glu19-Ser184), with a C-terminal 6-His tag

Product Description

Endotoxin Content:

Note: This product is for in vitro research use only and is not intended for use in humans or animals.