

## VPS54 Antibody

Catalog No: #42840



Package Size: #42840-1 50ul #42840-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

Product Name	VPS54 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antigen affinity purification.
Applications	IHC
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total VPS54 protein.
Immunogen Type	protein
Immunogen Description	Fusion protein of human VPS54
Target Name	VPS54
Other Names	WR; HCC8; SLP-8p; VPS54L; hVps54L; PPP1R164
Accession No.	Swiss-Prot#: Q9P1Q0Gene ID: 51542
Concentration	0.9mg/ml
Formulation	Rabbit IgG in pH7.4 PBS, 0.05% NaN <sub>3</sub> , 40% Glycerol.
Storage	Store at -20°C


## Application Details

Immunohistochemistry: 1:25-1:100

## Images



Immunohistochemical analysis of paraffin-embedded Human thyroid cancer tissue using #42840 at dilution 1/20.



Immunohistochemical analysis of paraffin-embedded Human liver cancer tissue using #42840 at dilution 1/20.

## Background

This gene encodes for a protein that in yeast forms part of a trimeric vacuolar-protein-sorting complex that is required for retrograde transport of proteins from prevacuoles to the late Golgi compartment. As in yeast, mammalian Vps54 proteins contain a coiled-coil region and dileucine motifs. Alternative splicing results in multiple transcript variants encoding different isoforms?

Note: This product is for in vitro research use only and is not intended for use in humans or animals.