

VPS37D Antibody

Catalog No: #42845



Package Size: #42845-1 50ul #42845-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

| | |
|-----------------------|---|
| Product Name | VPS37D Antibody |
| Host Species | Rabbit |
| Clonality | Polyclonal |
| Purification | Antigen affinity purification. |
| Applications | IHC |
| Species Reactivity | Hu |
| Specificity | The antibody detects endogenous levels of total VPS37D protein. |
| Immunogen Type | protein |
| Immunogen Description | Fusion protein of human VPS37D |
| Target Name | VPS37D |
| Other Names | WBSCR24 |
| Accession No. | Swiss-Prot#: Q86XT2Gene ID: 155382 |
| Concentration | 1.7mg/ml |
| Formulation | Rabbit IgG in pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol. |
| Storage | Store at -20°C |


Application Details

Immunohistochemistry: 1:25-1:100

Images



Immunohistochemical analysis of paraffin-embedded Human liver cancer tissue using #42845 at dilution 1/30.



Immunohistochemical analysis of paraffin-embedded Human thyroid cancer tissue using #42845 at dilution 1/30.

Background

Component of the ESCRT-I complex, a regulator of vesicular trafficking process. Required for the sorting of endocytic ubiquitinated cargos into multivesicular bodies. May be involved in cell growth and differentiation.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.