

## ZNF200 Antibody

Catalog No: #43209

Package Size: #43209-1 50ul #43209-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

Product Name	ZNF200 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antigen affinity purification.
Applications	IHC
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total ZNF200 protein.
Immunogen Description	Synthetic peptide of human ZNF200
Target Name	ZNF200
Other Names	MGC45293; Zinc finger protein 200; ZN200; ZNFMF
Accession No.	Swiss-Prot#: P98182Gene ID: 7752
Concentration	2mg/ml
Formulation	Rabbit IgG in pH7.4 PBS, 0.05% NaN <sub>3</sub> , 40% Glycerol.
Storage	Store at -20°C


## Application Details

Immunohistochemistry: 1:100-1:200

## Images



Immunohistochemical analysis of paraffin-embedded Human Lung tissue using #43209 at dilution 1/200.



Immunohistochemical analysis of paraffin-embedded Human Thyroid tissue using #43209 at dilution 1/200.

## Background

---

Could have a role in spermatogenesis.

---

Note: This product is for in vitro research use only and is not intended for use in humans or animals.