

## Recombinant Human Oncostatin M,OSM (N-6His)

Catalog No: #AP74837

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human Oncostatin M,OSM (N-6His)
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Accession No.	P13725
Calculated MW	24.44kDa
Target Sequence	MGSSHHHHHSSGLVPRGSHMAAIGSCSKEYRVLLGQLQKQTDLMQDTSRLLDPYIRIQGLDVPKLEHCRE RPGAFPSEETLRGLGRRGFLQTLNATLGCVLHRLADLEQRLPKAQLDERSGLNIEDLEKLQMARPNILGLRNNI YCMAQLLDNSDTAEPTKAGRGASQPPTPTPASDAFQRKLEGCRFLHGYHRFMHHSVGRVFSKWGESPNRSR R
Formulation	Lyophilized from a 0.2 µm filtered solution of 20mM PB, 150mM NaCl, pH 7.2.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks.Reconstituted protein solution can be stored at 4-7C for 2-7 days.Aliquots of reconstituted samples are stable at < -20C for 3 months.

## Application Details

tissue\_specificity:Homo sapiens (Human)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.