

Recombinant E. coli Ecotin (C-6His)

Catalog No: #AP74866

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant E. coli Ecotin (C-6His)
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Accession No.	P23827
Calculated MW	19.2kDa
Target Sequence	MKTILPAVLFAAFATTSAAWAESVQPLEKIAPYPQAEKGMKRQVIQLTPQEDESTLKVELLIGQTLEVDCNLHRL GGKLENKTLEGWGYDYVFDKVSSPVSTMMACPDGKKEKKFVTAYLGDAGMLRYNSKLPVVYTPDNVDVK YRVWKAEEKIDNAVVRLEHHHHHH
Formulation	Lyophilized from a 0.2 µm filtered solution of 20mM TrisHCl, 300mM NaCl, pH 8.5.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks.Reconstituted protein solution can be stored at 4-7C for 2-7 days.Aliquots of reconstituted samples are stable at < -20C for 3 months.

Application Details

tissue_specificity:Rattus norvegicus (Rat)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.