

Recombinant E. coli Methionine Aminopeptidase, MetAP, MAP (C-6His)

Catalog No: #AP74886

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant E. coli Methionine Aminopeptidase, MetAP, MAP (C-6His)
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Accession No.	P0AE18
Calculated MW	30.4kDa
Target Sequence	AISIKTPEDIEKMRVAGRLAAEVLEMIEPYVKPGVSTGELDRICNDYIVNEQHAVSACLGYHGYPKSVCISINEV CHGIPDDAKLLKDGDIVNIDVTVIKDGFGDTSKMFIVGKPTIMGERLCRITQESLYLALRMVKPGINLREIGAAI QKFVEAEGFSVVREYCGHGIGRGFHEEPQVLHYDSRETNVVLKPGMTFTIEPMVNAGKKEIRTMKDGWTVKT KDRSLSAQYEHTIVVTDNGCEILTLRKDDTIPAIISHDELEHHHHHH
Formulation	Supplied as a 0.2 µm filtered solution of 20mM TrisHCl, 100mM NaCl, 50% Glycerol, pH 8.0.
Storage	Store at < -20C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.

Application Details

tissue_specificity: Mus musculus (Mouse)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.