

# Recombinant E. coli G,U Mismatch-Specific DNA Glycosylase,Mug (C-6His)

Catalog No: #AP74887

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant E. coli G,U Mismatch-Specific DNA Glycosylase,Mug (C-6His)
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Accession No.	P0A9H1
Calculated MW	19.7kDa
Target Sequence	MVEDILAPGLRVVFCGINPGLSSAGTGFPFAHPANRFWKVIYQAGFTDRQLKPQEAQHLLDYRCGVTKLVDRP TVQANEVSKQELHAGGRKLIKIEDYQPQALAILGKQAYEQGFSQRGAQWGKQTLTIGSTQIWWLPNPSGLSR VSLEKLVEAYRELDQALVVRGRLEHHHHHH
Formulation	Supplied as a 0.2 µm filtered solution of 20mM TrisHCl, 2.5mM β-ME, 1mM PMSF, 50% Glycerol, pH 8.0.
Storage	Store at < -20C, stable for 6 months after receipt.Please minimize freeze-thaw cycles.

## Application Details

tissue\_specificity:Homo sapiens (Human)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.