

Recombinant Human Brain Natriuretic Peptide, BNP (His27-His134, N-6His)

Catalog No: #AP74889

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human Brain Natriuretic Peptide, BNP (His27-His134, N-6His)
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Accession No.	P16860
Calculated MW	14.2kDa
Target Sequence	MGSSHHHHHSSGLVPRGSHMHLGSPGSASDLETSGLQEQRNHLQGKLSLQVEQTSLEPLQESPRPTGV WKSREVATEGIRGHRKMLVLYTLRAPRSPKMQGSGCFGRKMDRISSSSGLGCKVLRH
Formulation	Supplied as a 0.2 µm filtered solution of 20mM PB, 150mM NaCl, 1mM EDTA, 1mM DTT, 20% Glycerol, pH 7.4.
Storage	Store at < -20C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.

Application Details

tissue_specificity: Mus musculus (Mouse)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.