

Recombinant Human Homeobox Protein OTX2 (C-6His)

Catalog No: #AP74890

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human Homeobox Protein OTX2 (C-6His)
Brief Description	Recombinant Protein
Host Species	E. coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Accession No.	P32243
Calculated MW	33.47kDa
Target Sequence	MMSYLKQPPYAVNGLSLTTSGMDLLHPSVGYPGPWASCPAATPRKQRRERTTFTRAQLDVLEALFAKTRYPD IFMREEVALKINLPESRVQVWFKNRRRAKCRQQQQQQNGGQNKVRPAKKTSPAREVSSESSTSGQFTPPS STSVPTIASSSAPVSIWSPASISPLSDPLSTSSSCMQRSYPMTYTQASGYSQGYAGSTSYFGGMDCGSYLTP MHHQLPGPGATLSPMGTNAVTSHLNQSPASLSTQGYGASSLGFNSTTDCLDYKDQTASWKLNFNADCLDYK DQTSSWKQVLEHHHHHH
Formulation	Lyophilized from a 0.2 µm filtered solution of 20mM PB, 150mM NaCl, pH 7.4.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7C for 2-7 days. Aliquots of reconstituted samples are stable at < -20C for 3 months.

Application Details

tissue_specificity: Tetronarce californica (Pacific electric ray) (Torpedo californica)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.