

Recombinant *S. cerevisiae* Translation Associated Element 1, TAE1, NTM1 (C-6His)

Catalog No: #AP74926

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant <i>S. cerevisiae</i> Translation Associated Element 1, TAE1, NTM1 (C-6His)
Brief Description	Recombinant Protein
Host Species	<i>E. coli</i>
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Accession No.	P38340
Calculated MW	27.13kDa
Target Sequence	MDVPADSHIKYEDAIDYWTDVDATVDGVLGGYGEGTVVPTMDVLGSNNFLRKLKSRMLPQENNVKYAVDIGA GIGRVSKTMLHKHAAKIDLVEPVKPFIEQMHVELAELKDKGQIGQIYEVGMQDWTPDAGKYWLIWCQWCVGH LPDAELVAFLKRCIVGLQPNGTIVVKENNTPTD DDFDETSSVTRSDAKFRQIFEEAGLKLASERQRGLPREL YPVRMYALKPMPNLEHHHHHH
Formulation	Supplied as a 0.2 µm filtered solution of 20mM PB, 150mM NaCl, pH 7.4.
Storage	Store at < -20C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.

Application Details

tissue_specificity: *Mus musculus* (Mouse)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.