

Recombinant Human S100 Calcium Binding Protein P,S100-P (N-6His)

Catalog No: #AP74990



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Recombinant Human S100 Calcium Binding Protein P,S100-P (N-6His)
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Accession No.	P25815
Calculated MW	12.56kDa
Target Sequence	MGSSHHHHHSSGLVPRGSHMTELETAMGMIIDVFSRYSGSEGSTQTLTKGELKVLMEKELPGFLQSGKDKD AVDKLLKDL DANGDAQVDFSEFIVFVAAITSACHKYFEKAGLK
Formulation	Supplied as a 0.2 µm filtered solution of 20mM Tris, 1mM DTT, 0.1M NaCl, 20% Glycerol, pH 8.0 .
Storage	Store at < -20C, stable for 6 months after receipt.Please minimize freeze-thaw cycles.

Application Details

tissue_specificity:Homo sapiens (Human)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.