

Recombinant E. coli N-Acetyl-D-Glucosamine Kinase,NAGK (C-6His)

Catalog No: #AP75303

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant E. coli N-Acetyl-D-Glucosamine Kinase,NAGK (C-6His)
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Accession No.	P75959
Calculated MW	34.5kDa
Target Sequence	MNHKVVHHHHHMYYGFDIGGKIALGVFDSGRQLQWEKRVPTPRDSYDAFLDAVCELVAEADQRFGCKGSV GIGIPGMPETEDGTLYAANVPAASGKPLRADLSARLDRDVRDNDANCFALSEAWDDEFTQYPLVMGLILGTG VGGGLIFNGKPITGKSYITGEFGHMRLPVDALTMMLDFPLRRCGCGQHGCNIENYLSGRGFAWLYQHYYHQP LQAPEIIALYDQGDEQARAHVERYLDLLAVCLGNILTVDPDLVIGGGLSNFPAITTLADRLPRHLLPVARVPR IERARHGDAGGMRGAAFLHLTD
Formulation	Supplied as a 0.2 µm filtered solution of 20mM TrisHCl, 100mM NaCl,2mM DTT,20% Glycerol, pH 8.0.
Storage	Store at < -20C, stable for 6 months after receipt.Please minimize freeze-thaw cycles.

Application Details

tissue_specificity:Homo sapiens (Human)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.