

# Recombinant E. coli 4'-Phosphopantetheinyl Transferase ACPS

Catalog No: #AP75918

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant E. coli 4'-Phosphopantetheinyl Transferase ACPS
Brief Description	Recombinant Protein
Host Species	E. coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Accession No.	P24224
Calculated MW	14.1kDa
Target Sequence	AILGLGTDIVEIARIEAVIARSGDRLARRVLSDNEWAIWKTHHQPVRFLAKRFVKEAAAKAFGTGIRNGLAFNQ FEVFNDELGKPRRLRWGEALKLAEKLGVANMHVTLADERHYACATVIIES
Formulation	Lyophilized from a 0.2 µm filtered solution of 20mM Tris, 150mM NaCl, pH 8.0.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7C for 2-7 days. Aliquots of reconstituted samples are stable at < -20C for 3 months.

## Application Details

tissue\_specificity:Homo sapiens (Human)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.