

Recombinant Human Endophilin-B1(SH3GLB1)

Catalog No: #AP76765

Package Size: #AP76765-1 20ug #AP76765-2 100ug #AP76765-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human Endophilin-B1(SH3GLB1)
Host Species	E.coli
Purification	Greater than 90% as determined by SDS-PAGE.
Immunogen Description	Expression Region:1-365aaSequence Info:Full Length
Target Name	SH3GLB1
Other Names	Bax-interacting factor 1
Accession No.	Q9Y371
Calculated MW	67.8 kDa
Tag Info	N-terminal GST-tagged
Target Sequence	MNIMDFNVKLAADAGTFLSRAVQFTEEKLGQAECTELDAHLENLLSKAECTKIWTEKIMKQTEVLLQPNNAR IEEFVYEKLDKAPSRINPELLGQYIMIDAGTEFGPGTAYGNALIKCGETQKRIGTADRELIQTSALNFLPLRN FIEGDYKTIAKERKLLQNKRLDLDAAKTRLKAKAAETRNSSEQLRITQSEFDRQAEITRLLLEGISSTHAHHLR CLNDFVEAQMTYYAQCYYMLDLQQLGSPSNYLSNNNQTSVTPVPSVLPNAIGSSAMASTSGLVITSPSNL SDLKECSGSRKARVLYDYDAANSTELSLLADEVITVFSVVGMDSDWLMGERGNQKGVPIYLELLN
Formulation	Tris-based buffer50% glycerol
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C,-80°C. The shelf life of lyophilized form is 12 months at -20°C,-80°C.Notes:Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Application Details

tissue_specificity:Homo sapiens (Human)

Product Description

1

Note: This product is for in vitro research use only and is not intended for use in humans or animals.