

## VWCE Antibody

Catalog No: #43431

Package Size: #43431-1 50ul #43431-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

Product Name	VWCE Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antigen affinity purification.
Applications	IHC
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total VWCE protein.
Immunogen Description	Synthetic peptide of human VWCE
Target Name	VWCE
Other Names	VWC1; URG11
Accession No.	Swiss-Prot#: Q96DN2 Gene ID: 220001
Concentration	1.8mg/ml
Formulation	Rabbit IgG in pH7.4 PBS, 0.05% NaN <sub>3</sub> , 40% Glycerol.
Storage	Store at -20°C


## Application Details

Immunohistochemistry: 1:20-1:100

## Images



Immunohistochemical analysis of paraffin-embedded Human esophagus cancer tissue using #43431 at dilution 1/30.



Immunohistochemical analysis of paraffin-embedded Human gastric cancer tissue using #43431 at dilution 1/30.

## Background

---

May be a regulatory element in the beta-catenin signaling pathway and a target for chemoprevention of hepatocellular carcinoma.

---

Note: This product is for in vitro research use only and is not intended for use in humans or animals.