

## ZNF473 Antibody

Catalog No: #43439

Package Size: #43439-1 50ul #43439-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

Product Name	ZNF473 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antigen affinity purification.
Applications	IHC
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total ZNF473 protein.
Immunogen Description	Synthetic peptide of human ZNF473
Target Name	ZNF473
Other Names	ZN473; HZFP100
Accession No.	Swiss-Prot#: Q8WTR7 Gene ID: 25888
Concentration	1.0mg/ml
Formulation	Rabbit IgG in pH7.4 PBS, 0.05% NaN <sub>3</sub> , 40% Glycerol.
Storage	Store at -20°C


## Application Details

Immunohistochemistry: 1:100-1:200

## Images



Immunohistochemical analysis of paraffin-embedded Human Liver cancer tissue using #43439 at dilution 1/200.



Immunohistochemical analysis of paraffin-embedded Human Lung cancer tissue using #43439 at dilution 1/200.

## Background

Involved in histone 3'-end pre-mRNA processing by associating with U7 snRNP and interacting with SLBP/pre-mRNA complex. Increases histone 3'-end pre-mRNA processing but has no effect on U7 snRNP levels, when overexpressed. Required for cell cycle progression from G1 to S phases.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.