

Recombinant Bovine Pregnancy-associated protein bPAP

Catalog No: #AP78015

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Bovine Pregnancy-associated protein bPAP
Brief Description	Recombinant Protein
Host Species	Baculovirus
Purification	Greater than 85% as determined by SDS-PAGE.
Immunogen Description	Expression Region:1-100aaSequence Info:Full Length
Target Name	N,A
Accession No.	P84291
Calculated MW	54.3 kDa
Tag Info	N-terminal MBP-tagged and C-terminal 6xHis-tagged
Target Sequence	DSELAGPRGARGPHGLSGPHGLSGLXGPXGYTGPIGMXGLTGLRRESEKVVLESKDGQELELVSSGSAQE ELELVSSGSAQVSFASYLGASQPLPSELW
Formulation	Tris-based buffer50% glycerol
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C,-80°C. The shelf life of lyophilized form is 12 months at -20°C,-80°C.Notes:Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Application Details

tissue_specificity:Hepatitis delta virus genotype I (isolate Italian) (HDV)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.