

## Glutamate--cysteine ligase catalytic subunit

Catalog No: #AP78990

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

|                       |  |
|-----------------------|--|
| Product Name          | Glutamate--cysteine ligase catalytic subunit   |
| Brief Description     | Recombinant Protein  |
| Host Species          | E.coli   |
| Purification          | Greater than 90% by SDS-PAGE   |
| Species Reactivity    | Human  |
| Immunogen Description | C-terminal   |
| Other Names           | GLCL, GLCLC, GCS heavy chain, Gamma-ECS, Gamma-glutamylcysteine synthetase   |
| Accession No.         | P48506 Gene name: GCLC   |
| Tag Info              | His  |
| Formulation           | 50mM NaH <sub>2</sub> PO <sub>4</sub> , 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (PH8.0)   |
| Storage               | Store at -20C. (Avoid repeated freezing and thawing.) Repeated freezing and thawing is not recommended.<br>Store working aliquots at 4°C for up to one week. |

Note: This product is for in vitro research use only and is not intended for use in humans or animals.