

Leucine-rich repeat-containing G-protein coupled receptor 4

Catalog No: #AP79726



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Leucine-rich repeat-containing G-protein coupled receptor 4
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	G-protein coupled receptor 48,GPR48
Accession No.	Uniprot ID:Q9BXB1
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.