

Inactive tyrosine-protein kinase transmembrane receptor ROR1

Catalog No: #AP79732



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Inactive tyrosine-protein kinase transmembrane receptor ROR1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Neurotrophic tyrosine kinase, receptor-related 1, NTRKR1
Accession No.	Uniprot ID:Q01973
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.