

Gamma-interferon-inducible lysosomal thiol reductase

Catalog No: #AP79744



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Gamma-interferon-inducible lysosomal thiol reductase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Gamma-interferon-inducible protein IP-30,Legumaturain,GILT,IP30
Accession No.	Uniprot ID:P13284
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.