

Phosphoglycerate kinase 1

Catalog No: #AP79748



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Phosphoglycerate kinase 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Cell migration-inducing gene 10 protein,Primer recognition protein 2,PRP 2,PGKA,MIG10,OK/SW-cl.110
Accession No.	Uniprot ID:P00558
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.