

Pulmonary surfactant-associated protein B

Catalog No: #AP79756



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Pulmonary surfactant-associated protein B
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	SP-B,6 kDa protein,Pulmonary surfactant-associated proteolipid SPL(Phe),SFTP3
Accession No.	Uniprot ID:P15781
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.