

# Programmed cell death 1 ligand 1

Catalog No: #AP79764



Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Programmed cell death 1 ligand 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	PD-L1,B7 homolog 1,B7-H1,CD274,B7H1,PDCD1L1,PDCD1LG1,PDL1
Accession No.	Uniprot ID:Q9NZQ7
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.