

RuvB-like 2

Catalog No: #AP79770



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	RuvB-like 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	48 kDa TATA box-binding protein-interacting protein,48 kDa TBP-interacting protein,51 kDa erythrocyte cytosolic protein,ECP-51,INO80 complex subunit J,Repressing pontin 52,Reptin 52,TIP49b,TIP60-associated protein 54-beta,TAP54-beta,INO80J,TIP48,TIP49B,CGI-46
Accession No.	Uniprot ID:Q9Y230
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.